	VitaCyte, LLC	Version: 3
	Product Insert	Date: 17 December 2021
Collagenase Gold Plus GMP Grade	Cat# 011-2000	

1. PRODUCT DESCRIPTION

Collagenase Gold Plus GMP Grade (Gold Plus) is an aseptically filled, lyophilized preparation of purified enriched *Clostridium histolyticum* collagenases. Gold Plus is sold as \approx 450 mg protein powder.

2. APPLICATION

Gold Plus is prepared with a minimum of 1500 FALGPA units per bottle¹. The required activity to recover cells will vary significantly depending on the tissue source and protocol used. Contact VitaCyte for technical guidance on how to evaluate the use of Gold Plus in specific applications. However, a concentration range of 0.5 – 1.5 mg/mL is a realistic starting point to evaluate cell recovery on many tissue types. Gold Plus contains minimal quantity of neutral protease activity and in the majority of cases will need to be supplemented with a neutral protease for successful cell recovery.

3. STORAGE & STABILITY

This product is stable for at least four years from date of manufacture if stored unopened as a lyophilized material at \leq 2-8°C. The product is shipped on Enviro ice packs to keep the product cold and minimize the potential for high temperature excursions during shipment.

4. PRODUCT USE

4.1. Enzyme Reconstitution


Reconstitute the lyophilized enzyme with at least 10 mL of water or buffer on ice for a minimum of 15 minutes to ensure complete dissolution of the enzyme. Occasionally invert the vial to aid in the dissolution process. The enzyme solution should not be vortexed or swirled excessively as enzyme denaturation may occur. Failure to allow the enzyme to completely rehydrate will affect the enzyme potency and could negatively impact the success of the tissue dissociation procedure. The enzyme is lyophilized in a buffer containing calcium so the initial reconstitution has sufficient calcium for enzyme stability. However, for optimal stability the final working buffer for tissue dissociation should have at least 0.1 mM Ca²⁺.

4.2. Digestion Solution Preparation

Once completely in solution, Gold Plus must be combined with a protease and diluted to the appropriate volume for use in a specific tissue dissociation procedure. To minimize loss of enzyme potency, the enzymes should be mixed just prior to beginning the digestion. At most, the mixture can be stored for two hours between 2°C and 8°C prior to use. This enzyme solution can be sterile filtered through 0.2 μ m cellulose acetate or PES filter membranes without compromising enzyme potency. Surfactant free cellulose acetate (SFCA) and PES filters from several major vendors were tested and no measurable loss of neutral protease activity was observed. The exact concentration of collagenase and protease is dependent on the specific application. Guidance for common cell targets is available at www.vitacyte.com.

5. TROUBLESHOOTING

5.1. Many factors contribute to the successful isolation of cells from tissue and inadvertent oversight to any of these conditions may drastically reduce the yield and viability of target cell population. While

	VitaCyte, LLC	Version: 3
	Product Insert	Date: 17 December 2021
Collagenase Gold Plus GMP Grade	Cat# 011-2000	

far from a complete list, the guidance below may help identify commonly encountered problems. Contact VitaCyte if this guidance does not help resolve specific issues.

5.2. Prolonged or Incomplete Digestion may be caused by:

- Loss of enzyme potency (activity)
- Incomplete enzyme rehydration during reconstitution
- Inappropriate enzyme dilution
- Presence of enzyme inhibitors
- Low incubation temperature

5.3. Low Yield and/or Cell Viability

- Prolonged organ/tissue warm ischemia time
- Aggressive mechanical disruption
- Extended incubation time
- Elevated incubation temperature
- Inappropriate enzyme dilution

6. ADDITIONAL INFORMATION

6.1. Intended Use & Regulatory


Gold Plus is for ex-vivo use only to recover cells from tissue. Guidance for use of reagents in clinical cell transplantation procedures is governed by local Institutional Review Boards and regional Health Authorities. Gold Plus is manufactured in accordance with the principles for clinical trial material outlined in Guidance Document Q7 Good Manufacturing Practice Guidance for Active Pharmaceutical Ingredients published by the FDA and ICH in September 2016. This product is labeled as ‘GMP Grade’ to indicate the quality system and manufacturing facility are consistent with requirements set forth in the above document. Guidance for the qualification and acceptance of reagents used in cell therapy applications can be found in the USP General Chapter <1043> Ancillary Materials for Cell, Gene, and Tissue-Engineered Products. Further details can be found on the VitaCyte Commitment to Quality document available upon request

6.2. Animal Origin

No bovine derived animal products are used in any step of manufacturing Gold Plus. Collagenase is purified from culture supernatants of *C. histolyticum* that contain porcine gelatin derived from European pigs. All other media components and the purification process are all entirely animal origin free.

6.3. Manufacturing Summary

Enzymes are purified from the culture supernatants results from the fermentation of native organisms. The purification processes use standard protein column chromatography and tangential flow filtration concentration and diafiltration techniques. The purification processes have been optimized to yield the highest purity attainable for each enzyme while minimizing undefined and contaminating protease activities. After characterization, the purified Gold Plus is sterile filtered in a

	VitaCyte, LLC	Version: 3
	Product Insert	Date: 17 December 2021
Collagenase Gold Plus GMP Grade	Cat# 011-2000	

qualified biosafety cabinet and aseptically dispensed into amber vials on activity units, lyophilized, sealed under vacuum then secured and labeled. The final lyophilized product is then further characterized to confirm each batch meets established specification ranges.

6.4. Activity Assessment

Each lot of product is characterized for collagenase activity using the FALGPA peptide substrate¹ and collagen degradation activity using fluorescein isothiocyanate conjugated to bovine Type I collagen fibrils². Endotoxin contamination is assessed using the Charles River Endosafe assay.

6.5. Additional Considerations

In addition to the quality of the dissociation enzymes, additional factors impact the outcome of success of cell isolations including: the quality of the organ/tissue and experience of the cell isolation team. The team needs to assess many variables that affect islet recovery. These include but are not limited to the characteristics of the donor, transport of the organ/tissue, the cell isolation procedure, and subsequent cell culture.

6.6. Resources & Support

Further details on manufacturing, quality control testing and use of products are available at www.vitacyte.com or technical support at 317-917-3457.

6.7. References

1. 1. Van Wart HE and Steinbrink DR. *Analytical Biochemistry* 113 (1981); 356-65.
2. 2. McCarthy RC et al. *Transplantation Proceedings* 40 (2008) 339-42.