



Pure science, defined.

SPECIFICATIONS
CIzyme™ BP Protease, Catalog # 003-1000
Revised June 2009

SPECIFICATIONS

Product Description	Purified neutral protease from <i>Bacillus polymyxa</i>
Lot #	
Appearance	White lyophilized cake in 10 mL amber bottle
Total Neutral Protease Activity ¹	> 800,000 NP Units
Endotoxin EU/mg protein ²	< 50 EU/mg

STORAGE: < -20 °C

EXPIRATION DATE: 2 Years from date of manufacture

¹ Voss EW, Workman CJ Mummert. Detection of protease activity using a fluorescence-enhancement globular substrate. *Biotechniques* 20:286-291 (1996)

² Lonza's Limulus Amebocyte Assay