



Pure science, defined.

PRODUCT SPECIFICATIONS
Cizyme™ Collagenase MA, Catalog # 001-2030
September 2009

SPECIFICATIONS

Product Description	Mixture of purified Class I and Class II Collagenase from <i>C. histolyticum</i>
Lot #	
Appearance	White lyophilized cake in amber bottle
Specific Collagen Degrading Activity (CDA) ¹	15,000 – 30,000 Units/mg
Total CDA Units per Container ¹	2.25 – 2.75 Million Units
Endotoxin EU/mg protein ²	< 10 EU/mg

STORAGE: < -20 °C

EXPIRATION DATE: 2 Years from date of manufacture

¹ McCarthy RC, et. al. Development and Characterization of a Collagen Degradation Assay to Assess Purified Collagenase Used in Islet Isolation. *Transplantation Proceedings* 40 (2008); 339-342

² Charles River Endosafe® PTS Endotoxin Assay